


Return Instructions:



1 Place vacutainers and cheek swab directly into specimen bag and seal. **DO NOT REMOVE ABSORBENT PAD.**

Place lab specimen bag into return mailer and seal.



3 Drop your return mailer bag off at a UPS Store **THE SAME DAY** and ask for a receipt with tracking number for proof of shipment and tracking purposes.



DO NOT place your return mailer bag inside a stand-alone drop box.



UPS Store



UPS Kiosk



Drug Testing With Dignity

Register Your Kit




ATTENTION: Please stop and scan our QR code on the back of your pouch to claim and register your Sharetek test. This is the mandatory first step to associate your results with desired name or alias. This step is **REQUIRED** to process and receive results from our lab. Thank you!

Collection Supplies

Kit Packet Includes:



Collection cup



Two vacutainer vials



Red tamper proof



Specimen bag



Cheek swab collection



Return mailer bag



Pouch bag

